- GB** **Instructions for use**
Targon® locking nail systems, in sterile packaging
- D** **Gebruuchsanweisung**
Targon®-Verriegelungsnagelsysteme, steril verpackt
- F** **Mode d'emploi**
Système de clous de verrouillage Targon®, sous conditionnement stérile
- E** **Instrucciones de manejo**
Sistemas de clavos de bloqueo Targon®, en envase estéril
- I** **Istruzioni per l'uso**
Sistemi di chiodi di bloccaggio Targon®, confezionati sterili
- P** **Instruções de utilização**
Sistemas de travamento Targon® em embalagem esterilizada
- NL** **Gebruiksaanwijzing**
Targon®-vergrendelingsnagelsystemen, steriel verpakt
- S** **Bruksanvisning**
Targon®-låsstiftsystem, sterilt förpackat
- RUS** **Инструкция по применению**
Системы фиксирующих гвоздей Targon®, стерильно упакованы
- CZ** **Návod k použití**
Blokovací systémy Targon®, sterilně balené
- PL** **Instrukcja użytkowania**
Systemy gwoździ ryglowanych Targon®, pakowanych sterylnie
- SK** **Návod na používanie**
Targon®-Systém uzatvárania klincov, sterilné balenie
- TR** **Kullanım Kılavuzu**
Targon® kilitlemeli iğne sistemi, steril paketinde

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 010481 03/12 V6 Änd.-Nr. 44403

Technical alterations reserved

Aesculap® Targon® locking nail systems, in sterile packaging

Intended use

The implant system is used for the intramedullary splinting, stabilization, and fixation of fractures in long tubular bones.

System configurations:

- Targon® PH/H
 - Proximal humerus nail
 - Humerus shaft nail
- Targon® PF, Targon® PFT
 - Proximal femur nail
- Targon® PF, Targon® TX
 - Tibia nail

Each system consists of locking nails, locking components (e.g., screws, sleeves, pins) and associated special implantation instruments.

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_T Titanium forged alloy Ti6Al4V acc. to ISO5832-3
- ISOTAN®_P Pure titanium acc. to ISO 5832-2
- PEEK-OPTIMA®

The titanium implants are coated with a colored oxide. Slight changes in coloration may occur, but do not affect the implant quality.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

PEEK-OPTIMA® is a registered trademark of Invivo, Ltd Lancashire FY5 4QD / UK.

Indications

Use Targon® PH/H for:

- Multi-fragment fractures of the humerus head, up to 4-fragment fracture (Targon® PH)
- Combined humeral head and shaft fractures (Targon® PH long)
- Stable or unstable humeral shaft fractures (Targon® H)
- Pathological fractures
- Pseudarthroses

Use Targon® PF/Targon® PFT for:

- Pertrochanteric, intertrochanteric and subtrochanteric femur fractures
- Above-mentioned fractures combined with ipsilateral femoral shaft fractures (Targon® PF long)
- Lateral femoral neck fractures combined with ipsilateral femoral shaft fractures (Targon® PF long)
- Peritrochanteric fractures (lateral femoral neck fractures extending into the trochanter zones)
- Pathological fractures
- Pseudarthroses

Use Targon® TX for:

- Open and closed tibial shaft fractures
- Repositioning osteotomies in the shaft region
- Reconstruction following tumor resection in the shaft region
- Osteotomies in the shaft region to lengthen or shorten the limb
- Pathological fractures
- Pseudarthroses
- Metaphyseal fractures of the tibia without articular involvement

The operating surgeon is responsible for establishing the above indications and indications beyond those, taking into account the specific clinical, biological, and biomechanical situation.

Contraindications

Do not use Targon® PH/H in the presence of:

- Fractures up to 7 cm proximal of the olecranon
- Crushed head cap of the humerus head (C3 fractures)
- Isolated avulsion fracture of a tuberosity

Do not use Targon® PF/Targon® PFT in the presence of:

- Medial and intermediary leg fractures

Do not use Targon® TX as the sole implant system in cases of:

- Fractures with articular involvement

Generally, do not use in cases of:

- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Anticipated overloading of the implant
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Side effects and interactions

- Dislocation, loosening and fracture of implant components
- Position change and loosening of fragments
- Delayed healing, or non-union and development of pseudarthrosis
- Primary and secondary infections
- Venous thrombosis, lung embolism, cardiac arrest
- Tissue reaction to implant materials
- Trauma to nerves, tendons, or vessels
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Compartment syndrome
- Pain in the region of the nail insertion point and in the region of locking components

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.

- Use implant components only in the nail holes intended for this purpose.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- Prior to the operation, the operating surgeon must familiarize him/herself with the Targon® implants and associated operating technique.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- The Targon® implants may only be implanted with the Aesculap implantation instruments intended for this purpose.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Mobilization, dynamic locking or use of compression involve the risk of migration of the nail and/or the locking component into the joint region or out of the bone. Therefore the migration paths of the implants and those of the fragments must be considered in the implantation and taken into account when selecting the implant length.
- Considerably higher strain due to increased leverage is to be expected if there are locking holes at the height or in the region of the fracture line. This can lead to implant failure. Accordingly, postoperative strain has to be reduced, or full load may be envisaged only after complete callus formation.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudarthrosis develops, or if the implant is overloaded for an excessively long period. One way to achieve this is through mobilization.
- Once the nail is locked, remove transmedullary bone and fragment correction screws again, or take into account, when determining the allowable patient load, that adverse lever effects can expose the nail to increased load.
- Always insert the Targon® implants in such a way that forces to be transmitted are low and can be transmitted by the bone at an early stage.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Always use cutting instruments, e.g., drill bits, guide pins, that are sharp.
- Do not apply strong force when implanting or explanting the implants. Should a problem arise, check the position and orientation of the implants, instruments, and fragments and analyze the origin of the problem. If necessary, repeat previous operating steps and check the instruments (e.g. for blocked drill helix).

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Increase the load on the fracture and implant in accordance with the progress of healing.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- At the time of follow-up exams, the scope of postoperative load and aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. Also, the dimensions of the implant should be taken into account.
- Recommended follow-up intervals
 - prior to the patient's discharge from hospital
 - 10-12 weeks post operationem
 - 6 months post operationem
 - 12 months post operationem
- Implanting the nail via a guide pin
 - Use an appropriate guide pin
 - Avoid jamming the implant or guide pin
 - Monitor the guide pin end by X-ray control
 - Remove the guide pin before locking the implant

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- ▶ Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- ▶ Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- ▶ Do not use implant components that are past their expiration date or whose packaging is damaged

Prior to resterilization, the implant components must be cleaned using the following validated reprocessing procedure:

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

Sterilization should always be favored over disinfection under all circumstances. Only if final sterilization is not possible should disinfection with a virucidal agent be preferred.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.aesculap-extra.net

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Validated reprocessing procedure

Mechanical alkaline cleaning and thermal disinfection

- ▶ Process the implant in its system storage device.
 - ▶ Place the implants on a tray that is suitable for cleaning (avoiding rinsing blind spots).
- Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - < 5 % anionic surfactant ■ 0.5 % working solution <ul style="list-style-type: none"> - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: B. Braun Helimatic Cleaner alkaline

For implant components that are to be resterilized:

WARNING

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- ▶ Handle the implants with new gloves only.
- ▶ Keep the implant system storage devices covered or closed.
- ▶ Process implant system storage devices separately from instrument trays.
- ▶ Clean implants must not be processed together with contaminated implants.
- ▶ Process the implant components individually and separately if no implant system storage devices are available, ensuring that the implant components are not damaged in the process.
- ▶ Mechanically clean and disinfect the implant components.
- ▶ Do not reuse surgically contaminated implants!

WARNING

Direct or indirect contamination may render implants unsuitable for resterilization!

- ▶ Do not reprocess implants that have been directly or indirectly contaminated with blood.

Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Immediately sort out damaged or inoperative products.

Packaging

- ▶ Place the product in its holder or on a suitable tray.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Make sure that the packaging will prevent a recontamination of the product whilst in storage.

Sterilization

- ▶ Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Application

WARNING

Risk of injury caused by incorrect operation of the product!

- ▶ Attend appropriate product training before using the product.
- ▶ For information about product training, please contact your national B. Braun/Aesculap agency.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- Treatment of a fracture with implants may not always be successful in completely restoring the original bone anatomy.
- Following the fracture treatment, the functionality of adjacent joints may be restricted.
- Pains may occur following the fracture treatment.
- The implants must not be exposed to overload due to excessive strains, hard physical work or sports activities. Overstraining entails the risk of implant loosening or breakage.
- The patient must be informed about the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The implants will not sustain the load and functionality of a healthy bone. The risks of transgressing these rules must be explained to the patient.
- If the implant loosens, fragments migrate, or pseudoarthrosis develops, a revision operation may be required.
- The patient must undergo regular medical follow-up examinations.
- The attending physician decides when to remove, in part or completely, the osteosynthesis implants through another surgical procedure.

WARNING

Complications due to incorrect use and/or incorrect assessment/treatment of the fracture!

- ▶ Make certain that nails and screws of a suitable length and diameter are used.
- ▶ Make certain that implants are aligned correctly.

- ▶ Check for correct mounting of the nail on the targeting device: Insert the drill bit through the drill sleeve and the targeting hole and guide it through the nail hole.

WARNING

Risk of nail seizure, nail bending, or bone rupture!

- ▶ Do not apply excessive force to instruments or implants.
- ▶ Apply only gentle hammer taps to insert Targon® TX and, if necessary, PF/PFT nails into the medullary cavity.
- ▶ Targon® To avoid bone rupture, do not hammer down PH/H nails under any circumstances.
- ▶ Check and, if necessary, correct the nail insertion point.
- ▶ If a nail proves difficult to insert, slightly drill open the medullary cavity.

The operating surgeon determines the following, depending upon the indications:

- whether implant locking is necessary
- if the implant is to be locked, whether it will be locked statically or dynamically.

Note

For locking an implant without a targeting attachment, Aesculap recommends free-hand locking supported by a C-bow and laser pilot beam, or with a radiolucent angle transmission!

WARNING

Avoid drill mistargeting, drill bit breakage, and guide pin breakage as well as premature implant failure due to crack formation as a result of damaged surfaces in the drill hole regions!

- ▶ Use a drill bit of the correct diameter.
- ▶ Level the drill insertion site with a flat counterbore or use a center punch.
- ▶ Only use sharp drill bits and guide pins, and drill carefully and without applying excessive force. Avoid direction deviations.
- ▶ During center punching, flat countersinking and drilling, do not apply force to the targeting device; do not bend or twist the targeting device.

Implant removal

- ▶ Prior to removing the nail, remove all screws.
- ▶ Remove the nail with quick, light taps, and with the aid of the impactor, extraction pin and slotted hammer.

Note

Implant removal can give rise to complications due to implants that are stuck, ingrown bone tissue, etc. Such complications can lead to damage to implants and/or instruments. For such cases we recommend using a set of special instruments, which can be ordered from Aesculap.

This instrument set should be at hand for every explanation. Other special instruments are available for dealing with broken implants. Note the instructions for use for such cases!

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

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