Vasco®

Instructions for Use

Glove for single use compliant with PPE regulation (EU) 2016/425, EN 420 and EN ISO 374

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Instructions for use are to be used in combination with the specific product related information on each product packaging. The gloves are sold as a bundled unit within the shipping carton. In case this bundled unit is dismantled and products are sold separately, the distributor must ensure that the instructions for use are accompanied with each separate unit.

The gloves are classified as Personal Protective Equipment (PPE) Category III according to PPE Regulation (EU) 2016/425 and have been shown to comply with this regulation through the applicable harmonised European standards. These gloves are designed to provide protection against specific chemicals tested, micro-organisms and particulate radioactive contamination (if applicable). The gloves meet the EN/ISO standards shown on each specific packaging. This PPE is single-use only and to be disposed after contamination.

EXPLANATION OF STANDARDS AND PICTOGRAMS

EN ISO 374-1	Permeation levels are based on breakthrough times (tested acc. EN 16523-1:2015) as follows:						
Type A / B / C	Permeation level acc. EN ISO 374-1:2016 +A1:2018	1	2	3	4	5	6
	Minimum breakthrough time in minutes	10	30	60	120	240	480
	Type A = chemical breakthrough time >30 minutes against at least 6 chemicals from the list Type B = chemical breakthrough time >30 minutes against at least 3 chemicals from the list Type C = chemical breakthrough time >10 minutes against at least 1 chemical from the list TEST CHEMICALS: A = Methanol / B = Acetone / C = Acetonitrile / D = Dichloromethane / E = Carbon disulphide / F = Toluene / G = Diethylamine / H = Tetrahydrofuran / I = Ethyl acetate / J = n-Heptane / K = Sodium hydroxide 40% / L = Sulphuric acid 96% / M = Nitric acid 65% / N = Acetic acid 99% / O = Ammonium hydroxide 25% / P = Hydrogen peroxide 30% / S = Hydrofluoric acid 40% / T = Formaldehyde 37%						
ABCDEFGHIJKLMNPST							
EN 374-4:2013	The degradation (in %) indicates the change in puncture resistance of the gloves after exposure to the respective challenge chemical.						

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections.

EN ISO 374-5: 2016	Tested for resistance to penetration according to EN 374-2:2014	EN 421:2010	Protection against particulate radioactive contamination.
	Tested for resistance to penetration by blood-borne pathogens according to EN ISO 374-5 / ASTM F1671 (virus resistance). Resistance to bacteria and fungi – pass Resistance to virus – pass The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.		These gloves do not protect against mechanical risks.
Virus			PPE is for single-use only and must not be reused.
 XXXX	XXXX = Identification number of notified Body responsible for the EU type examination and supervising ongoing conformity.		Before usage read instructions for use carefully.
		EN 420:2003+A1:2009	

PRECAUTIONS FOR USE

Always check the gloves for possible mechanical damage, e.g. holes or tears, before use. Do not use damaged gloves. Glove length is appropriate to the end use where the risk to the wrist area is minimal.

INGREDIENTS / HAZARDOUS INGREDIENTS

Some gloves might contain ingredients which are known to possibly cause skin irritations or allergic reactions with sensitised persons. Check warning information on specific packaging carefully. Formulation available on request.

TEMPERATURE LIMIT

STORAGE INSTRUCTIONS

Keep storage area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Protect gloves against ultraviolet light sources, sunlight, oxidizing agents and ozone. Store in original packaging according to the temperature limit, provided on the packaging.

DISPOSAL INSTRUCTIONS

Dispose of the gloves in accordance with the valid regulations for this material. Gloves contaminated with chemical substances must be disposed of in accordance with the regulations for the relevant chemicals.

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Title: PPER Info Sheet.pdf Initiator: Anke ? Fichtner

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Fichtner, Anke (fichande)
Title: OM-RA Manager RA
Date: Wednesday, 10 April 2019, 10:35 W. Europe Daylight Time
Meaning: Document signed as Author

UserName: Pausch, Gudrun (pausgude)
Title: OM-MD-DE08A - Director Global Marketing
Date: Wednesday, 10 April 2019, 13:33 W. Europe Daylight Time
Meaning: Approve Document

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tr KULLANIM KILAVUZU KKD Yönetmeliği (AB) 2016/425, EN 420 ve ISO 374'e uygun tek kullanımlık eldiven.